


Content:

- Safely reduce antibiotic exposure



You can also order free copies of additional material at info.pct@thermofisher.com

- Berlin consensus algorithms
- 5-minute abstracts of PCT key studies

For more content and free on-demand webinars please visit: thermoscientific.com/procalcitonin



Safely reduce antibiotic exposure

B·R·A·H·M·S Procalcitonin (PCT):
An effective tool for antibiotic stewardship




The challenge

Antibiotic resistance – an increasing threat to public health

Antibiotics (ABx) are a limited resource. At the current pace of injudicious use, all antibiotics will soon become ineffective. The WHO Global Action Plan on antimicrobial resistance, 2015, emphasizes that antimicrobial resistance is a crisis that must be managed with the utmost urgency.¹


The emergence and spread of antibiotic-resistant bacteria harm individuals and societies worldwide by causing:

- Prolonged illnesses
- Higher health care expenditures
- Greater risk of death



1/3 antibiotic prescriptions are unnecessary²

How does resistance to antibiotics develop?



A potential for change

B·R·A·H·M·S Procalcitonin (PCT) supports responsible use of antibiotics to prolong their effectiveness


Surviving Sepsis Campaign

International Guidelines for Management of Sepsis and Septic Shock, 2016³

“ We suggest that procalcitonin levels can be used to ...
... support **shortening the duration of antimicrobial therapy** in sepsis patients ...
... support the **discontinuation of empiric antibiotics** in patients who initially appeared to have sepsis, but subsequently have limited clinical evidence of infection.”

B·R·A·H·M·S
PCT

PCT-guidance of AB therapy has the potential to



Reduce initial prescription rates



Shorten treatment durations



Save overall treatment costs

Use of B·R·A·H·M·S PCT reduces antibiotic exposure

Strong evidence supports safe reduction of antibiotics using PCT-guided antibiotic stewardship protocols

- Reproducible, randomized clinical trials with more than 10,000 patients included
- Proven utility across diverse clinical settings: ICU, ED, Pediatrics, Neonatology, Surgery

Proven efficacy:
-16% to -74%
antibiotic exposure

No adverse impact on outcome






Figure 1 Relative reduction in AB exposure with PCT-guidance
 AB exposure in control group is normalized to 100, shown by a blue bar at the top. The gray bar depicts the relative exposure in PCT group and the orange bar shows the relative AB exposure reduction. All studies reported significant reduction in AB exposure.





When to start antibiotics?

B·R·A·H·M·S PCT enables rapid and reliable diagnosis of systemic bacterial infections¹⁹



PCT levels increase 3-6 hours after bacterial challenge and return to normal as the infection is resolved (Figure 2)^{19,20,21}

- ▶ High specificity and sensitivity for bacterial infection
- ▶ Indicator for disease severity and treatment response





Figure 2
Kinetics of PCT^{19,20}



How long to give antibiotics?

B·R·A·H·M·S PCT algorithms help tailor therapy to individual patient needs



PCT values should always be interpreted in context of the patient's clinical condition. Antibiotic treatment should be started/continued on suspicion of infection.

Daily monitoring of PCT course allows for customized ABx treatment duration, hence reduced ABx exposure

Ensure using the quality assay for SAFE clinical decision making

PCT cut-offs and clinical algorithms were established by use of the global reference standard Thermo Scientific™ B·R·A·H·M·S PCT™ sensitive KRYPTOR™ assay and are valid solely for all B·R·A·H·M·S PCT assays.



PCT-guidance for antibiotic therapy is a safe strategy

B·R·A·H·M·S PCT-guided reduction in antibiotic exposure could also reduce mortality rates




Figure 3 28-day mortality in PCT-guided group as compared to routine practice

Forest plot showing the comparison of PCT-guided algorithms vs. routine practice. The size of each square represents the proportion of information provided by each study. The vertical line depicts the point of “no difference” between the two groups, and the horizontal lines correspond to the 95% confidence intervals (CIs). Diamond represents the pooled odds ratio (OR) for all studies.

* 6-week follow-up

** 30-day follow-up

B·R·A·H·M·S PCT-guided antibiotic discontinuation Higher probability of survival

The lower mortality in PCT-guided patients may be attributed to

- Adequacy of antibiotics
- More timely recognition of alternative diagnoses
- Lower toxicity of antibiotics¹¹

- Randomized controlled interventional trial
- 1575 critically ill patients
- 15 centers





Figure 4 Probability of survival to day 365 in the PCT-guided group vs standard of care group¹¹

**Antibiotics are a double-edged sword.
Adequate dose helps, excess harms.**



Survival Benefit

Adults with LRTI symptoms




Patients in the ED

Is it bacterial infection?

As much as 75% of all antibiotic doses are prescribed for acute respiratory-tract infections, despite their mainly viral cause.⁵ PCT-guidance in such patients allows reduction in AB exposure without any adverse impact on outcome.⁴



-14%
initiation
of ABx⁴



-30%
adverse effects
from ABx⁴

Data from: Effect of Procalcitonin-Based Guidelines vs Standard Guidelines on Antibiotic Use in Lower Respiratory Tract Infections (ProHOSP)⁴

Largest prospective, multicentre, randomized controlled trial with PCT in LRTI patients presenting to EDs:


- 1359 LRTI patients, 6 centers
- PCT group (n=671), control group (n=688)

* % reduction related to non PCT-guided group


B·R·A·H·M·S PCT algorithm for LRTI patients



When to start ABx?



When to stop ABx?



PCT values should always be interpreted in context of the patient's clinical condition. Antibiotic treatment should be started/continued on suspicion of infection.

Community-acquired pneumonia (CAP)

Tailor the treatment duration in hospitalized patients





Figure 5 Meta-analysis data for 2027 patients hospitalized for CAP, total exposure of ABx in median days: PCT group = 6 days, control group = 10 days²⁴

Acute COPD exacerbations

Does every exacerbation require ABx?

- Significant sustained reduction in total antibiotic exposure for up to 6 months²⁵
- No decrease in mean time to next exacerbation²⁵
- No increase in lung function decline²⁵



If it is viral, antibiotics will not help. PCT can quickly identify patients who will benefit from antibiotic therapy.

Adults in Intensive Care Units



How to know the appropriateness of an empiric antibiotic?

Effective antibiotic treatment is reflected by declining PCT values,²⁶ consistent with its half-life time of about 20-24 hours.²⁰ Serial determinations of PCT can be used to monitor the course of infection in sepsis patients. Appropriate empiric antibiotic therapy was associated with a significant decline in PCT from day 2 to day 3 ($\Delta\text{PCT} \geq 30\%$).²⁶





Figure 6 Typical course of PCT serum level according to patient's response to antibiotic treatment (n=109)²⁷

B·R·A·H·M·S PCT algorithm for sepsis patients



Efficacy and safety in **critically ill patients**




Figure 7 Median duration of AB treatment in PCT-guided group = 5 days, in control group = 7 days¹¹



Data from: **The Stop Antibiotics on Procalcitonin Guidance Study¹¹**

- Largest prospective, multicentre, randomized, controlled, open-label intervention trial with PCT in critically ill patients
- Conducted in the Netherlands – a healthcare system with comparatively low use of ABx²⁸
- 1575 critically ill patients, 15 centers

Surgical ICU patients

Intra-abdominal infections are a common cause of infectious mortality in surgical ICUs. The duration of antibiotic treatment for their management is controversial.^{29,30}





Figure 8 Mean duration of AB treatment in PCT-guided group = 5 days, in control group = 10 days²⁹

**Excess of antibiotics is harmful.
PCT indicates the right time to stop.**

Children with LRTI symptoms




Children presenting to ED with **LRTI** – need for a targeted use of ABx

Antibiotics are overused in children and adolescents with LRTI.¹⁰
 PCT-guided treatment can markedly reduce ABx exposure in this patient group without any adverse impact on outcome.



Figure 9 Mean duration of AB treatment in PCT-guided group = 4.5 days, in control group = 6.3 days¹⁰

B·R·A·H·M·S PCT algorithm for children with LRTI symptoms



PCT values should always be interpreted in context of the patient's clinical condition. Antibiotic treatment should be started/continued on suspicion of infection.

Pediatric community-acquired pneumonia (CAP)

Pediatric CAP, in many cases, despite viral etiology is treated with antibiotics, leading to considerable over-use and increase in

- Risk of bacterial resistance
- Incidence of drug related adverse events
- Therapeutic costs⁹

PCT-guidance can help to avoid unnecessary antibiotics.



-14%
initiation
of ABx⁹

-85%
adverse effects
from ABx⁹




Figure 10 Duration of AB treatment in PCT-guided group = 5.37 days, in control group = 10.96 days⁹




Figure 11 Antibiotic exposure according to disease severity and treatment group⁹

Neonates with suspected early-onset sepsis




Early detection of neonatal sepsis

Avoid unnecessary ABx

Early diagnosis of neonatal sepsis is vital to improve the outcome. In the absence of reliable infection markers during the first hours of life, AB treatment in newborn infants with risk factors for infection is started early, exposing a considerable number of patients to unnecessary treatment.³¹

PCT-guidance has been shown to significantly reduce antibiotic treatment duration in such cases (Figure 12).¹⁷



B·R·A·H·M·S PCT cut-offs for neonates

To exclude maternal-fetal infection³²



PCT in cord blood <0.6 µg/L

Post-test probability of bacterial infection <0.001%

To reduce the duration of antibiotic treatment¹⁷



2 consecutive normal PCT values:

Stop antibiotic treatment for children born at term and prematurely

Consider PCT reference values according to children's age in hours (Figure 13)

In healthy neonates, plasma PCT concentrations increase gradually after birth, reaching peak values at about 24 hours of age and then decrease to normal values below 0.5 µg/L by 72 hours of age.¹⁷





Figure 13 Age adjusted PCT cut-off values in newborns¹⁷



-22.4 hours
antibiotic therapy





Figure 12 PCT-guided treatment of suspected neonatal sepsis¹⁷

PCT enables detection of neonatal sepsis from the first day of life.





The economic impact of PCT-guided treatment

has been studied through health economic modeling in various settings:

- Sepsis patients – ICU¹⁸
- Acute Respiratory Infections – inpatient, ICU, outpatient³³
- COPD exacerbation – inpatient³⁴

Treatment cost reductions ranging from 9% to 12% have been demonstrated across various countries.³⁵

The cost of testing for PCT is more than offset by downstream cost savings

“PCT helps me to prescribe antibiotics rationally and thus to save their power for future generations.”



PCT recognized by WHO as an aid to fight antibiotic resistance:

PCT included on the Essential In Vitro Diagnostics List of the World Health Organization to aid in antibiotic stewardship³⁷

References

- 1 WHO Global Action Plan Antimicrobial Resistance, 2015
- 2 Centers for Disease Control and Prevention (CDC), Press Release, May 3, 2016: <http://www.cdc.gov/media/releases/2016/p0503-unnecessary-prescriptions.html>
- 3 Rhodes et al., *Intensive Care Med* 2017; 43 (3): 304-377
- 4 Schuetz et al., *JAMA* 2009; 302: 1059-1066
- 5 Christ-Crain et al., *Lancet* 2004; 363: 600-607
- 6 Christ-Crain et al., *Am J Respir Crit Care Med* 2006; 174: 84-93
- 7 Long et al., *Crit Care* 2014; 18: 471
- 8 Briel et al., *Arch Intern Med* 2008; 168: 2000-7
- 9 Esposito et al., *Respir Med* 2011; 105: 1939-1945
- 10 Baer et al., *PLoS One* 2013; 8: e68419
- 11 de Jong et al., *Lancet Infect Dis* 2016; 3099: 1-9
- 12 Bouadma et al., *Lancet* 2010; 375: 463-474
- 13 Nobre et al., *Am J Respir Crit Care Med* 2008; 177: 498-505
- 14 Stolz et al., *Eur Respir J* 2009; 34: 1364-1375
- 15 Hochreiter et al., *Crit Care* 2009; 13: R83
- 16 Schroeder et al., *Langenbecks Arch Surg* 2009; 394: 221-226
- 17 Stocker et al., *Neonatology* 2010; 97: 165-174
- 18 Kip et al., *J Med Econ* 2015; 1-10
- 19 Harbarth et al., *Am J Respir Crit Care Med* 2001; 164: 396-402
- 20 Meisner M, *Procalcitonin – Biochemistry and Clinical Diagnosis*, ISBN 978-3-8374-1241-3, UNI-MED, Bremen 2010
- 21 Müller B et al., *Crit Care Med* 2000, 28 (4): 977-983
- 22 Svoboda et al., *Hepatogastroenterology* 2007; 54: 359-363
- 23 Maravić-Stojković et al., *Srp Arh Celok Lek* 2011; 139: 736-742
- 24 Schuetz et al., *Clin Infect Dis* 2012; 55(5): 651-62
- 25 Stolz et al., *Chest* 2007; 131: 9-19
- 26 Charles et al., *Crit Care* 2009; 13: R38
- 27 Stüber, 21st Int. Congress of Intensive Care and Emergency Medicine (ISICEM), Brussels 2001
- 28 Filippini et al., *Health Policy* 2006; 78(1): 77-92
- 29 Maseda et al., *J Crit Care* 2015; 30: 537-542
- 30 Sartelli et al., *Surg Infect (Larchmt)* 2015; 17: 9-12
- 31 Joram et al., *Arch Dis Child Fetal Neonatal Ed* 2006; 91: F65-F66
- 32 Joram et al., *Eur J Clin Microbiol Infect Dis* 2011; 30(8): 1005-1013
- 33 Schuetz et al., *Clin Chem Lab Med* 2015; 53(4): 583-592
- 34 Maas et al., Poster, ISPOR 18th Annual European Congress, Milan, 2015
- 35 Steuten, Cost effectiveness of biomarker testing in sepsis patients; ISICEM Brussels, March 15-19, 2015
- 36 Broyles, *Open Forum Infect Dis* 2017; 4(4): ofx213
- 37 2nd WHO Model List of Essential In Vitro Diagnostics (WHO/MVP/EMP/2019.05)

B·R·A·H·M·S PCT – Secured clinical decision making independent of platform



B·R·A·H·M·S PCT sensitive KRYPTOR™
B·R·A·H·M·S PCT direct



ALINITY i B·R·A·H·M·S PCT
ARCHITECT B·R·A·H·M·S PCT



VIDAS® B·R·A·H·M·S PCT



LIAISON® B·R·A·H·M·S PCT II GEN



Lumipulse® G B·R·A·H·M·S PCT



VITROS® B·R·A·H·M·S PCT



ELECSYS® B·R·A·H·M·S PCT



ADVIA Centaur® B·R·A·H·M·S PCT
Atellica® IM B·R·A·H·M·S PCT



Automated sensitive assay



Point-of-Care-Test

Not all intended uses and applications mentioned in this printing are registered for every product in every country.

Clinical Diagnostics

Thermo Fisher Scientific
B·R·A·H·M·S GmbH
Neuendorfstr. 25
16761 Hennigsdorf
Germany

+49 (0)3302 883 0
+49 (0)3302 883 100 fax
info.pct@thermofisher.com
www.thermoscientific.com/brahms

Find out more at [thermoscientific.com/procalcitonin](https://www.thermoscientific.com/procalcitonin)



Not all products are CE marked or have 510(k) clearance for sale in the U.S. Availability of products in each country depends on local regulatory marketing authorization status.

© 2020 Thermo Fisher Scientific Inc. All rights reserved. B·R·A·H·M·S PCT and all other trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. ADVIA Centaur and Atellica IM are registered and protected trademarks belonging to Siemens Healthcare Diagnostics. ADVIA Centaur B·R·A·H·M·S PCT and Atellica IM B·R·A·H·M·S PCT are products of Siemens Healthcare Diagnostics licensed from Thermo Fisher Scientific. ALINITY i and ARCHITECT B·R·A·H·M·S PCT are products of Abbott licensed from Thermo Fisher Scientific. Elecsys is a registered and protected trademark belonging to Roche or one of its subsidiaries. Elecsys B·R·A·H·M·S PCT is a product of Roche licensed from Thermo Fisher Scientific. LIAISON is a registered and protected trademark belonging to DiaSorin S.p.A. LIAISON B·R·A·H·M·S PCT II GEN is a product of DiaSorin S.p.A. licensed from Thermo Fisher Scientific. Lumipulse is a registered trademark of Fujirebio Inc. in Japan and in other countries. Lumipulse G B·R·A·H·M·S PCT is a product of Fujirebio Inc. licensed from Thermo Fisher Scientific. VIDAS is a registered trademark of bioMérieux S.A. or one of its subsidiaries. VIDAS B·R·A·H·M·S PCT is a product of bioMérieux licensed from Thermo Fisher Scientific. VITROS is a trademark of Ortho Clinical Diagnostics. VITROS B·R·A·H·M·S PCT is a product of Ortho Clinical Diagnostics licensed from Thermo Fisher Scientific. KRYPTOR is a trademark of Cisbio Bioassays, licensed for use by B·R·A·H·M·S GmbH, a part of Thermo Fisher Scientific. Patents: www.brahms.de/patents