



HDMI, WiFi, etc..



Controller

Use smartphone, tablet, laptop or desktop computer to change scenarios and control vital signs

Local Monitor

Connect your TV, display or projector to [Controller] for monitoring on the site

90 96 96 97 20

Remote Monitor

Use tablet, laptop or desktop computer to monitor remotely







Vital Signs Monitor Simulator for Emergency Care Training

We use EmerSim® for our maternal emergency life saving training. It allows me to prepare a scenario and quickly change preset vital signs values during the training. It gives trainees a real sense of emergency.

- Dr. Kayo Osada, Tsukuba Central Hospital



EXAMPLE 1: Training

Replicate real-life situations and scenarios for hands-on emergency care training.

EXAMPLE 2: Seminar

For lectures, connect EmerSim[®] to a projector using HDMI cable to display for a large audience.

EXAMPLE 3: Field Training

Set up a field scenario such as ambulance transportation practice, using a portable charger.

Key Features

- ✓ Connect with mobile and desktop devices on cloud.
- ✓ Customizable scenarios with various preset.
- ✓ The preset vital signs values can be easily changed with a single tap.
- ✓ Multiple sound effect options for realistic training enhancement.
- ✓ Real-time control, endless possibilities.
- ✓ No special equipment required.

EmerSim® Waveform Varieties

Electrocardiogram

Normal - NSR

1°block

2°block (Mobitz I)

2°block (Mobitz II)

3°block

Asystole

Atrial fibrillation

Atrial flutter

Junctional

LBBB

RBBB

ST elevation

SVT

Torsades de pointes

VF

Atrial flutter

VT

Arterial oxygen saturation

Normal Poor perfusion

